**Manufacturer**  
SHOFU INC.  
11 Kamitakamatsu-cho, Fukuine,  
Higashiyama-ku, Kyoto 605-0983, Japan  
www.shofu.com

**EC (REP)**  
SHOFU DENTAL GmbH  
Am Brill 17, 40878 Ratingen, Germany  
SHOFU DENTAL CORPORATION  
1225 Stone Drive, San Marcos,  
California 92078-4059, USA

SHOFU DENTAL ASIA-PACIFIC PTE.LTD.  
10 Science Park Road, #03-12,  
The Alpha, Science Park II, Singapore 117684

**SHOFU VINTAGE Art**

**EN** Instructions for use  
**DE** Gebrauchsanleitung  
**FR** Mode d'emploi  
**NL** Gebruiksaanwijzing  
**ES** Instrucciones de uso  
**IT** Istruzioni per l'uso  
**RU** Инструкция по применению  
**PL** Instrukcja obsługi  
**RO** Instrucțiuni de utilizare  
**SV** Bruksanvisning  
**EL** Οδηγίες χρήσης

02.2016-09  
70203

## EN Fluorescent Porcelain Stains

KEEP FOR FUTURE REFERENCE

VINTAGE Art are fluorescent porcelain stains for the internal and external modification of shades for high fusing porcelains.

### COLOR

Bl, O, W, B, Y, Gr, DR-Br, K, AS, BS, CS, G, Bl-G, CY, Br, B-Br, O-Br, F-AS, F-BS, F-CS, F-DS, MP, Mlv, P, V, DS, RS, RP, WR, P-O, F-RS, GP

### INDICATIONS

Shade adjustment of ceramic restorations

### PRECAUTIONS

- If any inflammation or other allergic reactions occur on operator, immediately discontinue use and seek medical advice.
- Avoid contact with skin and eyes. In case of accidental contact with skin, immediately blot with alcohol moistened cotton ball, and rise with plenty of water. In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
- Wear protective glasses etc while grinding and polishing this product to avoid damages to eyes.
- Use local dust extractor, dust protective mask, etc while grinding this product to avoid harmful influence of the dust on the human body.
- This product is intended for use by dental professionals only.

### DIRECTIONS FOR USE

- Cleaning of the restoration  
Clean the restoration thoroughly with water or alcohol etc to remove any dirt and other contaminants.
- Preparation  
Dispense the desired paste onto a glass plate or stain palette. The pastes can be mixed at any ratio to obtain the desired shade. The viscosity can be adjusted by mixing with VINTAGE Art STAIN LIQUID.
- Application  
Apply the paste with a brush.
- Firing

#### 4-1. For use of foundation or internal porcelain staining

Drying temp. <°C>	Drying time <min.>	Vacuum Start <°C>	Temp. Raising Speed <°C/min.>	Final temp <°C>	Vacuum End <°C>	Holding time <min.>
600	5-7	-	50-60	830-850	-	-

#### 4-2. For use of external staining

Drying temp. <°C>	Drying time <min.>	Vacuum Start <°C>	Temp. Raising Speed <°C/min.>	Final temp <°C>	Vacuum End <°C>	Holding time <min.>
600	5-7	-	50-60	850-870	-	-

#### 4-3. For characterizing of porcelain teeth

Drying temp. <°C>	Drying time <min.>	Vacuum Start <°C>	Temp. Raising Speed <°C/min.>	Final temp <°C>	Vacuum End <°C>	Holding time <min.>
600	5-7	-	50-60	1,000	-	2,0

### NOTE ON USE

- Tightly close the cap immediately after each use.
- Dispense necessary amount of paste and VINTAGE Art STAIN LIQUID. Do not put remaining material back into the container.
- Do not touch the material with bare hands.
- Use the dispensed material immediately after dispensing.
- Optimal firing conditions vary due to the different designs and operating voltages of the porcelain furnaces. It is essential to carry out test firings before firing actual restorations.
- Do not use any liquid other than the VINTAGE Art STAIN LIQUID.
- Follow the instructions for use of each dental material or instrument to be used in conjunction with this product.
- Avoid any source of ignition since this product is flammable.
- Do not mix with any other products.
- Do not use this product for any purpose other than specifically outlined in INDICATIONS in these instructions for use.

### STORAGE

Store at room temperature(1-30 °C/34-86 °F). Keep away from direct sunlight and any source of ignition.

**CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.**

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