

Glaslonomer Cement *CX-Plus* CAPSULE

For Predictable Permanent Cementation



GlasIonomer Cement CX-Plus CAPSULE

Simple, Hassle-free and Reliable ...

GlasIonomer CX - Plus Capsule

A time tested luting cement for –

Cementation of crown & bridge work, inlays and orthodontic appliances

Features

- Durable chemical bond to both tooth structure and metal
- Low film thickness of just 15 microns, optimal thixotropicity and ideal flow characteristics ensure a complete marginal seal
- High resistance to moisture contamination at the initial stage
- Radiopaque with high fluoride release
- 3 year shelf life

Physical Properties

Minimum Deliverable Amount (ml)		0.14
Working Time (23° C) From start of mixing		3' 10"
Setting Time (37° C) From start of mixing		4' 10"
Tensile Bond Strength (MPa)	To Enamel	5.8
	To Dentin	4.4
Film Thickness (μ)		15
Acid Erosion (mm)		0.02



Availability

- Packed in a box of 50's (0.53gm / capsule) [PN 3277]



SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605 - 0983, Japan (www.shofu.co.jp)

SHOFU DENTAL ASIA-PACIFIC PTE. LTD. 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)

SHOFU DENTAL TRADING (SHANGHAI) CO. LTD. No.645 Jiye Road, Sheshan Industry Park Songjiang Shanghai 201602, China (www.shofu.com.cn)

SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)

SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany (www.shofu.de)

SHOFU UK. Riverside House, River Lawn Road, Tonbridge, Kent TN9 1EP, United Kingdom (www.shofu.co.uk)