

GlasIonomer FX-II CAPSULE

Enhanced Direct Restorative



GlasIonomer FX-II CAPSULE



Fast, Easy and Predictable with every single mix

GlasIonomer FX-II Capsules

An enhanced general purpose direct restorative for

- Minimal intervention (MI) treatment
- Pediatric and geriatric restoration
- Intermediate restorative for heavy stress bearing Class I & II restorations
- Restorations of Class III & V cavities including abfraction defects
- Final restorative for non-stress bearing Class I and II cavities in adult dentition
- For core build-up and base / liner

Clinical Case 1



Clinical Case 2



Photographs courtesy of Kenji Murai, DDS, Ph.D

Features

- Durable mechanical properties
- A strong chemical bond to dentin and enamel with continuous fluoride release and recharge
- Biocompatible and radiopaque with low solubility
- Fluorescence similar to natural tooth structure
- 3 year shelf life



Physical Properties

Working Time (23 °C) [from start of mixing]	1' 40"
Setting Time (37 °C) [from start of mixing]	2' 30"

Availability

- Packed in a box of 50's
(0.55gm / capsule)
- Shades:
A2 [PN 3270], A3 [PN 3271],
A3.5 [PN 3272], B2 [PN 3273]



SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605 - 0983, Japan (www.shofu.co.jp)
 SHOFU DENTAL ASIA-PACIFIC PTE. LTD. 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)
 SHOFU DENTAL TRADING (SHANGHAI) CO. LTD. No.645 Jiye Road, Sheshan Industry Park Songjiang Shanghai 201602, China (www.shofu.com.cn)
 SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)
 SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany (www.shofu.de)
 SHOFU UK. Riverside House, River Lawn Road, Tonbridge, Kent TN9 1EP, United Kingdom (www.shofu.co.uk)