

# ENDURA ZERO

Flat Plane 0° Posteriors



**ENDURA ZERO** is a line extension of the 3-layered opalescent composite resin. Endura Posterio Teeth which offers, trusted quality with excellent toughness, strength and durability. The unique spillway and edge design of the 0° posterior teeth, increase masticatory efficiency and provides easy alignment for monoplane occlusion.



**Ideal for advanced ridge resorption & deviated mandibular position**

**Facilitates greater range of lateral movement**

**More comfort with minimal stress on supporting oral tissues**

## BENEFITS



### IMPROVED MASTICATORY FUNCTION

Straight edge contact during intercuspation with wider buccolingual table



### PREVENTS FOOD RETENTION

Unique spillway design in buccolingual and mesiodistal direction



### AESTHETICS & COMFORT

Natural flow of cervical line with smooth transition from anterior to posterior teeth

## MOULD CHART

**Shades: A3, A3.5**

Measurement [mm]	UPPER	Width of 4 posteriors Length of first premolar (minus collar)
		LOWER

**M28**

28.0

8.0

30.0

9.0



**M30**

30.0

8.5

32.5

9.5



**M32**

32.0


9.0

35.0

10.0




## ORDER INFORMATION



### ENDURA ZERO


Packed in a set of 8's  
(12 sets per box)

## RELATED PRODUCTS



### ENDURA ANTERIO

Packed in a set of 6's  
(16 sets per box)



### ENDURA POSTERIO

Packed in a set of 8's  
(12 sets per box)



**SHOFU DENTAL ASIA-PACIFIC PTE. LTD.**

1 Science Park Road, #01-01 / 02 The Capricorn, Science Park II,  
Singapore 117528 Tel: (65) 6377 2722  
eMail: mailbx@shofu.com.sg www.shofu.com.sg



ZRO20250103