



CLICK / SCAN



FOR MORE INFO

HIGH DENSITY HIGH STRENGTH AESTHETIC PRESS CERAMIC

A high-density, high-strength lithium disilicate ceramic designed with HDIL (High Density Interlocking) technology for highly aesthetic restorations with exceptional strength.

- **Remarkable durability** with a flexural strength of 525 MPa and excellent acid resistance
- **Stable shades** and antagonist friendly with precise margins
- Produces **minimal reaction layers** for fast finishing and a more efficient workflow
- **43 shades** in four levels of translucency to meet diverse aesthetic needs



Available in:

- Individual shades available as refills - 5 ingots, 3g each
- Trial Kit [REF P0870] 5 ingots, 3g each HT-BL4, HT-A1, HT-A2, LT-A1 and LT-A2

CLINICAL APPLICATIONS

- Veneers
- Inlays/Onlays
- Crowns/Partial Crowns
- 3 units Anterior bridge

SHADE SYSTEMS

Shade	Inlay/ Onlay	Veneer	Partial Crown	Crown	Bridge
E	●	●	●	-	-
HT	●	●	●	●	●
LT	●	●	●	●	●
HO	-	-	-	●	●

FEATURES

- ✓ **High Strength**
 - Exceptional flexural strength (525 MPa) ensures greater longevity.
- ✓ **Aesthetic Versatility**
 - Offers individualized aesthetics with Vintage LD layering porcelain and Vintage Art Universal stain
- ✓ **Colour Stability**
 - Retains shade and translucency even after multiple firings.
- ✓ **Biocompatibility & Intra-oral Stability**
 - Antagonist-friendly with excellent abrasion resistance.
 - High acid resistance with low solubility.

RELATED PRODUCTS

- Vintage Art Universal Stains
- Vintage LD Porcelain
- Ceravety Press & Cast
- Zilmaster Finishing & Polishing Kit



Official Partner

